

Application Instructions for Allflex Electronic Identification Ear Tags

Apply Allflex Electronic ID Tags with the Red Universal Total Tagger.

Use the **red** blunt pin and **remove** the black insert from the base of the jaws for all EID Tags.



1. To load, depress spring clip and insert the female EID tag. Ensure that the raised portion of the tag, which encloses the transponder chip, is placed in the open portion of the jaws.
2. Slip the male tag completely onto the blunt applicator pin. Squeeze the jaws together lightly to ensure the male shaft is in line with the female.
3. Dip the jaws of the applicator holding the tag into an antiseptic or disinfectant solution.
Suggestions:
Iodine Solution
Nolvasan
Vetericyn

Recommended EID Tag Placement

Application site must be free of foreign debris prior to placement of tags on the animal. Review application instructions prior to tagging.

IMPORTANT: Caution, "Free Air Space" is critical for proper healing and retention. Inspect placement after tagging to ensure there is sufficient space between ear and EID tag. It is recommended to pull tag pieces apart for proper healing.



For CATTLE—Left Ear Only

1. The EID tag should be placed in the middle of the ear between the two cartilage ribs close to the head.
2. The female (EID) portion of the tag should be on the inside of the ear.



For SWINE—Left Ear Preferred

1. Place the EID tag in the middle or thicker part of the ear.
2. The male tag should be on the outside of the ear and the female (EID) should be placed on the inside of the ear. Tag retention is improved when male tag lays flush against the back side of the ear.

IMPORTANT: Tag application may vary based on breed type and ear size.



For SHEEP & MEAT GOATS—Left Ear Preferred

1. Place the EID tag in between the cartilage ribs of the ear, near the first third of the ear (closest to the head).
2. The female EID portion of the tag should be placed on the inside of the ear.

IMPORTANT: Disinfect tags prior to tagging and spray for 3-5 days post-tagging to reduce the opportunity for infection.

4-H Beef Tagging Guidelines and Tips for EID Tags

Instructional videos for appropriate Allflex EID application may be found here:

<https://www.allflex.global/na/support/>

Applicators – for the EID tags, you **MUST** use the red **Allflex Universal Total Tagger** with the red blunt pin. *You must remove the black plastic insert from the base of the jaws.*

Placement of the EID tag in the applicator - The beef tags have a bump on the tag that is caused by the coil of wire. There is an arrow (line) on the other side of the tag. Place this to the top of the applicator when loading the tag. * See the pictures on the reverse side.

Application of the tags – As you can see in the instructions, the beef tags **MUST** be put in the **LEFT EAR** with the EID button facing to the front (inside the ear). They should be placed close to the head (1/4th to 1/3rd of the ear length) out from the base, in the middle of the ear vertically.

DO NOT put the tags in an existing hole in the ear. Doing so will greatly increase the chance of loss.

*See the pictures on the reverse side. After application, the EID tag should spin freely while the 4-H tag is held. If not, the male stem was not pushed far enough into the EID tag.

Use disinfectant to dip the tags and applicator in. It is recommended that you use a product such as Nolvasan, Vetericyn, Tek-Trol, or an iodine solution. They are effective, non-irritation products and can be purchased from the local vet or vet supply store. Using a disinfectant will help the tag site heal and provide lubrication for easier tag application. Using a product that irritates the ear (alcohol, bleach, etc.) will cause the animal to do more rubbing of the ear and increase the chance of losing the tag.

Keeping things warm – It does not matter when we schedule the steer tagging, it will be cold. Cold tags are harder to put together. Here are some ideas to keep them warm. You can keep your disinfectant in heated bucket or put the bucket of warm water inside cooler. The tags will go together easier if kept warm. If you have electricity, you can use a space heater to keep the tags in front of. If not, a cooler with plastic bottles of hot water in it can be used to keep the tags in.

Keep the matched EID tags and 4-H tags together – The tags are packaged with the EID and visual 4-H numbers matched and in a database available from Allflex. If for some reason the tags should get separated, you can use the sheet of labels to match them back up. Allflex has been working on their website over the past few years. If you are unable to download your tag data and need it, contact Lexie.

Using the labels – Each bag of tags has 2 sets of labels that have the EID tag, a barcode of that 15-digit number, and the visual 4-H number on it. You can use one set of labels for county tagging if you wish. You should keep one set of labels in a file to be used later for any state nominations. *You should also make a copy or take a picture of the labels and put them in a safe place, in the event your labels should get damaged or misplaced.*

Weigh-In Card – A template of the county weigh-in card can be downloaded from the KSU YLP website.

If the steer already has an EID tag – It is unlawful to remove an 840 tag from any species of livestock, as it is their official (ADT) identification. An animal may also only have one EID tag. The Kansas 4-H EID tags are official 840 EID tags. Please contact Lexie if you need to tag a 4-H/FFA project with an existing EID tag. We have a process for those animals being state nominated. **DO NOT** tag an animal with an additional 840 EID tag, if they already have one – it can legally only have one.



Questions – Can be directed to Lexie Hayes, KSU Youth Livestock Program Coordinator:

adhayes@ksu.edu or 785-532-1264

Place Line at Top



Place Bump at Top



Kansas 4-H Swine, Sheep, & Meat Goat Tagging Guidelines and Tips for EID Tags

We have experienced a good retention rate on the 4-H EID tags. We realize that sometimes “stuff happens”, but we can often trace a tag replacement issue to the application procedure. While it is impossible to have 100% tag retention, it can be increased by following these instructions and tips.

Instructional videos for appropriate Allflex EID application may be found here:

<https://www.allflex.global/na/support/>

Applicators – for the EID tags, you **MUST** use the red **Allflex Universal Total Tagger** with the red blunt pin. The black plastic insert needs to be removed for application, just as it was for the beef tags. If you still have small livestock tags that were distributed from 2013-2016, please contact Lexie. The tagger configuration and instructions are slightly different than they are for the tags we’ve used for the last few years.

Application of the tags – As you can see in the instructions, the tags are to be put in the ear with the EID button facing to the front (inside the ear). It is preferred that they be placed in the left ear (as you look at the animal from the rear). They should be placed close to the head (1/3 to 1/2 of the ear length out from the base), in the middle of the ear vertically. **DO NOT** put the tags in an existing hole in the ear. Doing so will greatly increase the chance of loss. Pulling on both sides of the tag to increase air flow around the tag and turning the button will verify correct application and ensure there is appropriate space between the tag and the animal’s ear.

Use disinfectant to dip the tags and applicator in. It is recommended that you use a product such as Nolvasan, Vetericyn, or an iodine solution. They are effective, non-irritating products that can be purchased from the local vet or vet supply store. Using a disinfectant will help the tag site heal and it will provide lubrication for easier tag application. **It will also reduce the opportunity for infection and is imperative to use disinfectant when tagging meat goats and sheep!** Using a product that irritates the ear (alcohol, bleach, etc.) will cause the animal to do more rubbing of the ear and increase the chance of losing the tag.

After Application – Grasp both parts of the tag and pull on them. This will insure that the stem is completely in and it will create room for air flow around the tag to promote faster healing.

Keep the matched EID button and 4-H visual tags together – The tags are packaged with the EID and visual 4-H numbers matched and in a database available from Allflex USA. If for some reason the tags should get separated, you can use the sheet of labels to match them back up. Allflex has been working on their website over the past few years. If you are unable to download your tag data and need it, contact Lexie.

Using the labels – Each bag of tags has 2 sets of labels that have the EID tag, a barcode of that 15-digit number, and the visual 4-H number on it. You can use one set of labels for county tagging if you wish. You should keep one set of labels in a file to be used later for any state nominations. **Please also make a copy of the labels and put them in a safe place, in the event your labels should get misplaced.**

Weigh-In Card - A template of a county weigh-in card can be downloaded from the KSU YLP website.

If the animal already has a EID 4-H tag – It is unlawful to remove an 840 tag from any specie of livestock, as it is their official (ADT) identification. An animal may also only have one EID. The Kansas 4-H EID tags are official 840 EID tags. Please contact Lexie if you need to tag a 4-H/FFA project with an existing EID tag. Any non-840 tags may be removed and replaced with the Kansas 4-H EID tag.

Questions – Can be directed to Lexie Hayes, KSU Youth Livestock Program Coordinator:

adhayes@ksu.edu or (785)532-1264